



COMPOSITE & CERAMIC POLISHERS - SINGLE USE

# DISPOSABLE KENDA NOBILIS

## 1-step diamond polishing system


For single patient use.


One polisher only for cost effective finishing and polishing of composite restorations and ceramics. The flexibility of the diamond polisher allows a gentle to teeth clean finishing and polishing of composite restorations and ceramics.




Perfect smooth surfaces will be achieved. Disposable polisher saves sterilisation cost and time.

Recommendation:

Disposable KENDA UNICUS 1-Step Diamond Polisher for ultimate high shine polishing (page K-8).



Order REF	PS 0309	PS 0311	PS 0312
Size Ø x L mm	10x1	6x7.5	5x10
ISO No.	808 204... 303 513 100	030 513 060	243 513 050
Grit	fine		
Color	violet		
	5000 RPM		

-  Box of: 6, 12, 25, 100
-  Pressure: 0,3 - 1 N
- Shank: RA, medical plastic
- Medical device: class IIa
-  disposable