



COMPOSITE & CERAMIC FINAL POLISHERS - SINGLE USE

DISPOSABLE KENDA UNICUS 

1-step diamond polishing system

For single patient use.

Ultra fine diamond polisher for final polishing of all composites and ceramics.


This polisher can be used on composites as well as ceramics without changing the polisher.

The benefits are time saving as well as cost saving.




Disposable polisher saves sterilisation cost and time.

Recommendation:

Disposable KENDA NOBILIS 1-Step Diamond Polisher for finishing and polishing (page K-6).



Order REF	PS 0359	PS 0361	PS0362
Size Ø x L mm	10x1	6x7.5	5x10
ISO No.	808 204...	303 492 100	030 492 060
Grit	very-fine		
Color	orange		
	5000 RPM		

-  Box of: 6, 12, 25, 100
-  Pressure: 0,3 - 1 N
- Shank: RA, medical plastic
- Medical device: class IIa
-  disposable